

**Product Name :**  
RTD Studies in Packed Bed**Product Code :**  
CHINAELABC0150001**Description :**

RTD Studies in Packed Bed

**Technical Specification :**

Rota meter is provided to measure the flow.

For understanding the R.T.D. characteristics, a special arrangement to inject tracer into the lower end of reactor (for liquid and gas both) is provided.

The setup consists of a glass column packed with Rasching Rings.

Liquid is fed to the reactor through liquid distributor, fitted at the bottom of the column.

The flow rate can be adjusted by operating the needle valve provided on Rota meter.

Samples can be taken periodically from the top outlet of reactor Experiments.

To plot RTD curve for Plug Flow Reactor.

To determine the Dispersion No Utilities Required Water Supply & Drain Instruments, Laboratory Glassware and Chemicals required for analysis as per the system adopted.

